



You decide when to fill and when to go

The iFill system makes it easy and convenient for you to fill portable oxygen cylinders right in your home to give you greater freedom to come and go on your terms.

- The iFill system gives you the flexibility to fill a wide range of cylinders from the M4 and M6, popular for traveling, to back-up E cylinders.
- Select a cylinder size that best suits your out-of-home experience. For example, lightweight M6 cylinders provide up to 4 hours of portable oxygen, giving you the option to leave home without an oxygen cart. An M6 cylinder, and others, can be comfortably carried by shoulder strap or around the waist.
- All cylinders can be refilled as often as you need.

Specifications (continued)

Oxygen Duration

NOTE: All ambulatory ranges are calculated assuming a breath rate of 20 BPM in PulseDose® mode.

Flow Rate:	Use Times (Shown In Hours)							
	1	1.5	2	2.5	3	4	5	6
M4								
PulseDose	5.7	3.8	2.9	2.3	1.9	1.4	1.1	.9
Continuous Flow	1.9	1.3	.9	.7	.6	.5	.4	.3
M6								
PulseDose	8.3	5.5	4.1	3.3	2.8	2.1	1.7	1.4
Continuous Flow	2.7	1.8	1.4	1.1	.9	.7	.6	.4
ML6								
PulseDose	8.6	5.7	4.3	3.4	2.9	2.1	1.7	1.4
Continuous Flow	2.8	1.9	1.4	1.1	.9	.7	.6	.5
C								
PulseDose	12.1	8.1	6.1	4.9	4.0	3.0	2.4	2.0
Continuous Flow	4.0	2.7	2.0	1.6	1.3	1.0	.8	.7
D								
PulseDose	21.0	14.0	10.5	8.4	7.0	5.2	4.2	3.5
Continuous Flow	6.9	4.6	3.5	2.8	2.3	1.7	1.4	1.2
E								
PulseDose	34.4	23.0	17.2	13.8	11.5	8.6	6.9	5.8
Continuous Flow	11.4	7.6	5.7	4.6	3.8	2.8	2.3	1.9

This chart is intended to be used only as a guide.



DeVilbiss Healthcare LLC
100 DeVilbiss Drive
Somerset, PA 15501

800.338.1988 • FAX 800.552.7379
www.DeVilbissHealthcare.com

©2015 DeVilbiss Healthcare LLC 11.15 All Rights Reserved.
DeVilbiss®, PulseDose® and DeVilbiss iFill® are registered trademarks of DeVilbiss Healthcare. LT-2279 Rev D

What you want is what you get

independence.
innovation.
iFill.



Independence on Your Terms

Home oxygen therapy is intended to enhance your quality of life, and encourage independence and activity. The iFill Personal Oxygen Station was created just for you, giving you greater independence, *on your terms*. No more waiting for oxygen deliveries or missing important functions; the iFill system lets you decide when you want to leave your home.





The DeVilbiss iFill® puts you in control

The DeVilbiss iFill Personal Oxygen Station represents the latest in portable oxygen convenience, providing unmatched flexibility and ease of use. The iFill unit extracts purified oxygen from room air using standard oxygen concentrator technology.*

- Existing oxygen patient? Don't worry – the iFill system can be used with any manufacturer's concentrator ... even yours.
- Connectors on the cylinder and iFill unit simply click into place. Control panel lights indicate when the cylinder is filling and when the cylinder is full.
- An M6 iFill cylinder can be filled in just a little over an hour – not the two hours it can take with other in-home oxygen filling units. Shorter fill times save energy costs.

* The unit is designed to fill oxygen cylinders only and cannot be used for direct oxygen inspiration.

- The iFill system can be placed and used virtually anywhere in the home – minimizing the noise and disruption to daily life.

Choice of cylinders

DeVilbiss PD1000A PulseDose® Oxygen Conserving Cylinder is widely recognized as the gold standard of conserving technology, offering true 3.5:1 savings. The PD1000A conserving device is designed to minimize discomfort associated with oxygen therapy. Patients report decreased drying and irritation of the nose.

- Quickly delivers a 16.5 ml bolus of oxygen (per liter setting) deep into your lungs. The volume of oxygen is consistent, regardless of your breathing rate, eliminating the need to adjust settings during exertion.
- Intuitive click-stop prescription control knob lets you easily choose among eight common prescription settings.
- Rugged, impact-resistant GE LEXAN® case.
- Single lumen cannula.
- Visual indicators for low battery (flashing red LED), insufficient battery (solid red LED) and normal pulsing (flashing green LED).
- Two standard AA alkaline or NiMH batteries provide approximately 50 days of operation (2 lpm prescription setting; 20 bpm, 4 hours/day). Continuous flow mode requires no batteries.


Standard Regulator

- Pneumatic continuous flow regulator
- Provides .5–6 lpm continuous flow
- Single lumen cannula
- Brass fittings

LEXAN is a registered trademark of General Electric Company.

iFill 870 Post Valve Cylinders

- Can be used with most conserving devices.



Specifications

DeVilbiss iFill Personal Oxygen Station

Output	93% O ₂ ± 3%
Width	12.25 inches
Height	28.5 inches
Depth	22.5 inches
Weight	66 lbs.
Warranty	3 years or 3000 hours – whichever comes first

Cylinder Fill Times


NOTE: All filling times are approximate and may vary depending on environmental conditions, such as altitude. Typical cylinder fill times to 2,000 ± 100 psig are as follows:

M4	60 minutes
M6	75 minutes
ML6	90 minutes
C	130 minutes
D	215 minutes
E	350 minutes

PD1000A Series Cylinder Weights

M4	3.5 lbs.
M6	4.1 lbs.
ML6	4.8 lbs.
C	5.6 lbs.
D	7.2 lbs.
E	9.8 lbs.

Specifications continued on back cover.



Get out and live

Research findings confirm that active patients on Continuous Oxygen Therapy (COT) have a significantly higher survival rate than inactive patients on Nocturnal Oxygen Therapy (NOT).

Source: Petty TL, Bliss PL, Ambulatory Oxygen Therapy, Exercise, and Survival with Advanced Chronic Obstructive Pulmonary Disease (The Nocturnal Oxygen Therapy Trial Revisited) Respiratory Care 2000; Vol 45(2) :208.

NOTE: All specifications are subject to change without notice.